




Maxiofacial Surgery Instruments

A Bard-Parker handle with a No. 15 blade is the usual scalpel used to make incisions on the gums where a flap can be raised to gain access to tooth.



Retractors

Tongue, cheek and flap retractors are needed to provide adequate access to the tooth and are legion in number.

Needle Holders and Sutures

Sutures are used to replace and secure the gum flap in its former position, to limit the size of the wound and therefore restrict subsequent contamination with food debris, and to assist in controlling bleeding. The usual suture size for intra-oral work is 4/0 and the material may be absorbable or non-absorbable Scissors
Dissecting scissors and suture-cutting scissors are usually used.



Artery Forceps

These are designed to hold soft tissue without damaging it. College tweezers are used to lift up sutures prior to removing them.

Periosteal Elevators

A number are available and are similarly used like a normal elevator but periosteal elevators have two heads – one at each end. One end is sharp-tipped to expand the tooth socket while the other has a blunt end to remove granulation tissues or raise a gum flap.



Chisels and Mallet

Chisels are more appropriate in the operating theatre and are particularly useful for disto-angular third molars and upper wisdom teeth extraction.

Bone Curettes

The Mitchells trimmer is probably the most valuable instrument in the category and is used to remove granulation tissue in the tooth socket.

Bone Ronguer and Bone File

Used to smoothen any sharp bone edges after tooth extraction is done.





- Роторасширители
- Ecarteurs - Ouvre-bouches
- Vestibulumabhalter / Mundspanner
- Роторозширювачі
- Separadores - Abre-bocas
- الكاشطات - الكمامات الفم



AMi-947
Mouth Gag Molt
11 cm for (Children)



AMi-948
Mouth Gag Molt
14 cm for (Adults)



Cheek Retractor
AMi-949 140 mm
AMi-950 180 mm



Cheek Retractor
AMi-951 140 mm



AMi-952
Wangen Holder,
vestibule
15 cm, with or without
lip protection,



AMi-953
Wangen Holder,
vestibule 15 cm, with
or without lip protection,



AMi-954
Wangen Holder
Wangen holder with dual
function, 17 cm,



AMi-955
140 mm

AMi-956
150 mm

Cheek Retractor



- Роторасширители
- Ecarteurs - Ouvre-bouches
- Vestibulumabhalter / Mundspanner
- Роторозширювачі
- Separadores - Abre-bocas
- الكاشطات - الكمامات الفم



AMi-957
Fig. 1
weitlander
10.5cm
13.5cm



AMi-958
Fig. 2
weitlander
10.5cm
13.5cm



AMi-959
Stainless Steel
Double
Ended Cheek
Retractors



AMi-960
Stainless Steel
Single
Ended Cheek
Retractors



Weitlander
AMi-961 11cm
AMi-962 13cm
AMi-963 15cm



AMi-967
Denhard
12 cm



AMi-968
Roser-König
16 cm



AMi-969
Roser-König
19 cm



Jennings
AMi-964 11cm
AMi-965 13cm
AMi-966 15cm





Wound, Bendable
Retractor, individually or as saving kit!

- Ранорасширители
- Ecarteurs
- Wundhaken

- Ранорозширювачі
- Separadores
- الكاشطات



AMi-970
Cheek Retractor
Long



AMi-971
Cheek Retractor
Long



AMi-972
Cheek Retractor
Long



Retractor with
Ergo Grip
AMi-973 24 cm,
ergonomic handle,



AMi-974
Cheek Retractors,
Retractor Ivory
Width adjustable
& fixable, 20 cm,



Cheek Retractors,
tongue depressors
AMi-975 19 cm



Cheek Retractors,
Minnesota
AMi-976 16 cm



Cheek Retractors,
AMi-977 20 cm
easy to bend
by hand,



Cheek Retractors,
retractor and tongue
pushers BISHOP with
comfortable margin,
AMi-978 16 cm,



Retractors Cheek
retractor combination
Combination,
AMi-979 25 cm,

RETRACTORS

ASIA CARE + TECH



- Ранорасширители
- Ecarteurs
- Wundhaken
- Ранорозширювачі
- Separadores
- الكاشطات



AMi-982
Retractors, Unique
Length 175 mm,
angled 45°



AMi-983
Retractors, Unique
Length 175 mm,
angled 90°



AMi-984
Retractors
16cm



AMi-985
Retractors
19cm



AMi-986
Retractors, Mini-Langenbeck
16 cm, hooks 22 mm,



AMi-987
Retractors, Cheek
retractors combination
21 cm, 2



AMi-988
Mini Retractor
about 8 cm,



AMi-989
Retractors
17 cm, hooks 22 mm





- Ранорасширители
- Ecarteurs
- Wundhaken

- Ранорозширювачі
- Separadores
- الكاشحات



30 x 16 mm

AMi-990
Retractors, surgical-flat
surgically-flat, 21 cm, hook
30x16 mm,



AMi-991
Retractors, surgical-flat
surgically-flat, 21 cm, hook
42x11 mm



60 x 12 mm

AMi-992
Retractors, surgical-flat
surgically-flat, 21 cm, hook
63x12 mm,



Retractors, surgical-flat
surgically-flat, 21 cm, hook
AMi-993 30x16 mm,
AMi-994 42x11 mm
AMi-995 63x12 mm,



Retractors, Langenbeck
21 cm, hook
AMi-996 30x11 mm
AMi-997 30x16 mm,



AMi-998
Retractors
21 cm,



AMi-999
Farabeuf
(120 X 12mm
Set of 2 pices)



AMi-1000
Retractor Seldin



AMi-1001
FELDMAN
RETRACTOR
5-1/2", 14cm



GUM SCISSORS

ASIA CARE + TECH



- Ножницы десневые
- Ciseaux a gencives
- Zahnfleischscheren
- Ножиці для ясен
- Tijeras para encias
- مقص صمغ



AMI-1002
Iris
115mm
Straight



AMI-1003
Iris
115mm
Curved



AMI-1004
Goldman-Fox
130mm
Saw edge, straight



AMI-1005
Goldman-Fox
130mm
Saw edge, Curved



AMI-1006
La Grange
115mm
Saw edge, double Curved



Mayo
AMI-1007 140mm
AMI-1008 170mm
AMI-1009 230mm
Straight



Mayo
AMI-1010 140mm
AMI-1011 170mm
AMI-1012 230mm
Curved



Metzenbaum
AMI-1013 145mm
AMI-1014 180mm
AMI-1015 200mm
Straight



Metzenbaum
AMI-1016 145mm
AMI-1017 180mm
AMI-1018 200mm
Curved



AMI-1019
Kelly
160mm
Straight



AMI-1020
Kelly
160mm
Curved



AMI-1021
Iris
115mm
Straight



AMI-1022
Iris
115mm
Curved



AMI-1023
La Grange
115mm
Double Curved



AMI-1024
Goldman-Fox
130mm, Straight



AMI-1025
Goldman-Fox
130mm, Curved



AMI-1026
Kelly
160mm, Curved



AMI-1027
Iris
115mm, Straight



AMI-1028
Iris
115mm, Curved



AMI-1029
Micro Iris
90mm, Curved





- Ножницы десневые
- Ciseaux a gencivas
- Zahnfleischscheren

- Ножиці для ясен
- Tijeras para encias
- مقص صمغ



AMI-1030
Iris
115mm
Angled



AMI-1031
Iris
115mm
Curved



AMI-1032
Goldman-Fox
130mm
Saw edge,
straight



AMI-1033
Goldman-Fox
130mm
Saw edge,
curved



AMI-1034
Goldman-Fox
130mm
Saw edge,
double curved



AMI-1035
La Grange
115mm
Saw edge,
double curved



AMI-1036
Neumann
120mm
Saw edge,
double curved



AMI-1037
Kilner (Ragnell)
120mm
Straight



AMI-1038
Kilner (Ragnell)
120mm
Curved



AMI-1039
Joseph
140mm
Straight



AMI-1040
Joseph
140mm
Curved



AMI-1041
Kelly
160mm
Straight



AMI-1042
Kelly
160mm
Curved




AMI-1043
Locklin
160mm
Saw edge,
angular



AMI-1044
Locklin
120mm
160mm
Saw edge,
angular



AMI-1045
Dean
170mm
Saw edge,
angular



AMI-1046
Spencer
90mm



AMI-1047
Spencer
130mm



AMI-1048
Spencer
115mm



AMI-1049
Noyes
110mm, Straight



AMI-1050
Noyes
110mm
Curved



AMI-1051
Noyes
110mm
Angular



AMI-1052
Castroviejo
150mm
Curved

NEEDLE HOLDERS

ASIA CARE + TECH



- Иглодержатели
- Porte-aiguilles
- Nadelhalter

- Голкотримачи
- Porta-agujas
- اصحاب ابرة



AMi-1053
130mm
Straight



AMi-1054 Mayo
130mm
AMi-1055 Mayo
160mm
AMi-1056 Mayo
180mm
Straight Curved



AMi-1057 Mayo
140mm
AMi-1058 Mayo
200mm
Straight Curved



AMi-1059 Metzenbaum
145mm
AMi-1060 Metzenbaum
200mm
Straight Curved



AMi-1061
Beebee
110mm
Saw edge,
straight




AMi-1062
Beebee
110mm
Saw edge,
curved



AMi-1063
Beebee
110mm
Saw edge, straight



AMi-1064
Beebee
110mm
Saw edge, curved



AMi-1065
Universal
120mm
Saw edge



AMi-1066
Universal
120mm



AMi-1067 Crile
130mm
AMi-1068 Crile
150mm
AMi-1069 Crile
180mm



AMi-1070
Halsey
130mm



AMi-1071
Baumgartner
130mm



AMi-1072
Derf
120mm



AMi-1073
Delicate
150mm



AMi-1074 Mayo-Hegar
130mm
AMi-1075 Mayo-Hegar
160mm
AMi-1076 Mayo-Hegar
180mm
AMi-1077 Mayo-Hegar
200mm



AMi-1078
Crile-Wood
150mm
straight



AMi-1079
Crile-Wood
150mm
curved





- Ножницы десневые
- Ciseaux a gencives
- Zahnfleischscheren
- Ножци для ясен
- Tijeras para encias
- مقص صمغ



AMi-04-1080
Goldman Fox Scissors
12.5cm Curved



AMi-04-1081
Iris Scissors
11.5cm Curved
Serrated For Slip-Free Cutting



AMi-04-1082
LaGrange Scissors
11.5cm Double Curved
Serrated For Slip-Free Cutting

Double Curved



AMi-04-1083
Castroviejo Micro Scissors
12cm Curved




AMi-04-1084
Super Cut Goldman Fox
13.2cm Curved
Serrated For Slip-Free Cutting



AMi-04-1085
Super Cut Iris
11cm Straight
Serrated For Slip-Free Cutting



AMi-04-1086
Super Cut Iris
11cm Curved
Serrated For Slip-Free Cutting



AMi-04-1087
Super Cut LaGrange
11.5cm
Double Curved
Serrated For Slip-Free Cutting

Double Curved



- Ножницы десневые
- Ciseaux a gencivas
- Zahnfleischscheren
- Ножиці для ясен
- Tijeras para encias
- مقص صمغ

T.C. / CERAMIC COATED SCISSORS



Mayo Dissecting Scissors



AMi-04-1088
14.5cm
AMi-04-1089
17cm
AMi-04-1090
23cm
B/B, Str.



AMi-04-1091
14.5cm
AMi-04-1092
17cm
AMi-04-1093
23cm
B/B, Cvd.



Metzenbaum Standard Scissors



AMi-04-1094
14.5cm
AMi-04-1095
18cm
AMi-04-1096
20cm
AMi-04-1097
23cm
AMi-04-1098
26cm
AMi-04-1099
28cm
AMi-04-1100
30cm
B/B, Str.



AMi-04-1101
14.5cm
AMi-04-1102
18cm
AMi-04-1103
20cm
AMi-04-1104
23cm
AMi-04-1105
26cm
AMi-04-1106
28cm
AMi-04-1107
30cm
B/B, Cvd.



Locklin Dissecting Scissors, S/S

AMi-04-1108
16cm
1 Blade
Serrated



Locklin Dissecting Scissors, S/S

AMi-04-1109
16cm
1 Blade
Serrated



Dean Dissecting Scissors, S/S

AMi-04-1110
17cm
1 Blade
Serrated



- Иглодержатели
- Porte-aiguilles
- Nadelhalter

- Голкотримачи
- Porta-agujas
- اصحاب ابرة



Castroviejo
AMi-1080 130mm
AMi-1081 150mm
AMi-1082 180mm



AMi-1083
 Halsey
 130mm



AMi-1084
 Baumgartner
 130mm



AMi-1085
 Derf
 120mm



AMi-1086
 Delicate
 150mm



AMi-1087
 Adson
 170mm



Mayo-Hegar
AMi-1088 130mm
AMi-1089 160mm
AMi-1090 180mm
AMi-1091 200mm



AMi-1092
 Crile-Wood
 150mm



Castroviejo
AMi-1093 str 140mm cvd
AMi-1094 str 180mm cvd



Olsen-Hegar Needle
 Holder With Scissors
AMi-1095 140mm



Gillies Needle
 Holder With Scissors
AMi-1096 160mm



Mathieu
AMi-1097 140mm
AMi-1098 170mm



Castroviejo
AMi-1099 str 140mm cvd
AMi-1100 str 180mm cvd



Olsen-Hegar Needle
 Holder With Scissors
AMi-1101 140mm



Gillies Needle
 Holder With Scissors
AMi-1102 160mm



Hartmann (Baby
 Mosquito)
AMi-1103 90mm straight




Halstead Mosquito
AMi-1105 125mm
 Straight

AMi-1106
 Curved

AMi-1104
 125mm
 Curved



- Зажимы артериальные, зажимы для салфеток
- Pincés hemostatiques, Pincés a compreses
- Arterienklemmen, Tuychklemmen
- Затискачі артеріальні, затискачі для серветок
- Pinzas hemostaticas, Pincas para compresas
- الملقط هيموستاتيک المشابك منشفة



Halstead Mosquito 1x2
AMi-1107 125mm
straight



Kelly
AMi-1109 140mm
straight



Rochester-Pean
AMi-1111 140mm
AMi-1112 160mm
AMi-1113 180mm
AMi-1114 200mm
straight



Rochester-Kocher 1x2
AMi-1119 140mm
AMi-1120 160mm
AMi-1121 180mm
AMi-1122 200mm
straight



AMi-1127
Kelly-Rankin
160mm
straight



Lorna
AMi-1129 100mm
AMi-1130 130mm



AMi-1131
Allis-Baby 3:4
150mm



AMi-1132
Allis 4:5
200mm



AMi-1133
Allis 6:7
200mm



Backhaus
AMi-1134 140mm
Not-perforating and
not-damaging



AMi-1135
Backhaus
130mm



AMi-1136
Backhaus
85mm



AMi-1137
Jones
90mm



Gross Maier
Dressing Forceps
AMi-1138 200mm straight



Foerster
Dressing Forceps
AMi-1140 200mm



Young
Tongue Forceps
AMi-1141 170mm



Collins
Tongue Forceps
AMi-1142 170mm





- Инструменты для системы коффердам
- Pincas porte-crampons a digue
- Kofferdam Klammerzangen
- Инструменти для системи коффердам
- Pinzas de engrapar para dique de caucho
- السد المطاط الملقط المشابك




AMi-1143

On request the Content of the rubber dam kit can be changed with other instruments.

Content:

- 01 Anisworth
- 02 Brewer
- 03 Clamps Stand
- 04 Denture Plate
- 05 Rubber Dam Frame
- 06 Rubber Dam Box Lower Case
- 07 Rubber Dam Box Upper Case
- 08 Rubber Dam Clamp Fig. 212
- 09 Rubber Dam Clamp Fig. 211
- 10 Rubber Dam Clamp Fig. 8
- 11 Rubber Dam Clamp Fig. 8A
- 12 Rubber Dam Clamp Fig. 14
- 13 Rubber Dam Clamp Fig. 14A
- 14 Rubber Dam Clamp Fig. 0
- 15 Rubber Dam Clamp Fig. 00
- 16 Rubber Dam Clamp Fig. 2
- 17 Rubber Dam Clamp Fig. 2A
- 18 Rubber Dam Clamp Fig. 7




AMI1-1143
Small
Rubber Dam
Frame



AMI2-1143
Medium
Rubber Dam
Frame



AMI3-1143
Large
Rubber Dam
Frame



AMi-1145
Brewer



AMi1-1145
Stockes



AMi-1147
Palmer



AMi-1146
Ivory light weight



AMi-1144
Anisworth



AMi1-1144
Anisworth



Punch Plier



AMi2-1144
Anisworth



AMi3-1144
Ivory



AMi-1148
Seaparating
Elastic Plier



RUBBER DAM CLAMPS

ASIA CARE + TECH



- Кламмера для коффердама
- Crampons pour digue
- Kofferdam Klammern
- Кламмера для коффердама
- Grapas para dique de caucho
- المشابك السد المطاط





- Инструменты для удаления коронок
- Arrache-couronnes
- Teleskopkronenzangen

- Инструменты для видалення коронок
- Arranca-corona
- مزيلات التاج



AMI-1210
Crown Spreader



AMI-1211
Mershon
14 cm
Band Pusher



AMI-1212 Manual Crown Removers
Fig. 0 Handle with 3 tips
Only 3 tips



AMI2-1212



AMI3-1212



AMI4-1212



AMI6-1212



AMI7-1212



AMI8-1212


Manual Crown Removers
Handle with 3 tips
Only 3 tips



AMI-1213
Manual Crown Remover



AMI-1214
Complete with 3 tips
and tightening keys



semi-automatic crown
& bridge remover



AMI-1215



AMI-1216
Manual Crown Remover



- Инструменты для удаления коронок
- Arrache-couronnes
- Teleskopkronenzangen

- Инструменти для видалення коронок
- Arranca-corona
- مزيلات التاج



AMi-04-1217
160 mm

Crown Plier with plastic inserts Straight
For ceramical and plastic crowns and bridges, with
plastic inserts (sterilization up to 180° C).



AMi-04-1218
145 mm

Crown Plier
with plastic inserts
Angled



AMi-04-1219
140 mm

Universal plier for the removal
of ceramic crowns and bridges



AMi-04-1220
140 mm

Universal plier for the removal
of ceramic crowns and bridges



AMi-04-1221
Furrer
Copper Ring
Removing Plier.
150 mm
Curved



AMi-04-1222 125 mm
AMi-04-1223 150 mm
Böhm
Pliers for telescope
crowns.
Diamond points coated.



AMi-04-1224
Nerve-Canal Pliers
straight



AMi-04-1225
Nerve-Canal Pliers
Curved



- Матрицедержатели
- Porte Matrices
- Matrizenspanner

- Матрицеутримувачи
- Porta Matrices
- مصفوفة الخدم



AMI-1217
Fig. 5
Ivory 8
5 cm



AMi-1218
Fig. 6
Ivory 8
6 cm



AMI-1219
Fig. 7
Ivory 8
7 cm



AMI-1220
Fig. 1
Tofflemire-Universal



AMI87-1220



AMI1-1220



AMI2-1220



AMI-1221
Fig. 2
Tofflemire-Junior



AMI-1222
Fig. 3
Tofflemire-Senior



AMI-1223
Narrow



AMI-1224
WiAMi

Siqveland Retainers

Amalgam Well



AMI-1225
Stainless Steel
Amalgam Well



AMI-1226
Ivory 12



AMI-1227
Ivory 1



AMI-1228
Ivory



AMI-1229
Elliot





- Инструменты для амальгамы
- Instruments pour amalgame
- Amalgamis\instrumente

- Инструменти для амальгами
- Instrumentos para amalgama
- الصوك المنعم

- Ручки для скальпелей
- Manches de bistouris
- Skalpellgriffe

- Ручки для скальпелів
- Mangos de bisturi
- مشرط مقابض



Ø 2mm

AMi-1230
Fig. 10
165mm
English pattern



Ø 2mm

AMi-1231
Fig. 1
185mm
German pattern



Ø 2mm

AMi-1232
Fig. 1
150mm
Swiss pattern



Ø 2mm

AMi-1233
Fig. 1
150mm



Ø 2mm

AMi-1234
Fig. 2
150mm, 45°



Ø 2mm

AMi-1235
Fig. 3
140mm, 90°



AMi-1236
Fig.3
5" (12.7cm)



AMi-1237
Fig.3G
5" (12.7cm)



AMi-1238
Fig.4
5 1/2" (13.7cm)



AMi-1239
6 Way Adjustable Angled Scalpel Handle



Locks At 90°, 60°, 30°,
Straight, 120°, 150°, 180°,
Angles



- Ручки для скальпелей
- Manches de bistouris
- Skalpellgriffe
- Ручки для скальпелів
- Mangos de bisturi
- مشرط مقابض



Scalpel Handles

Stainless steel scalpel handle



AMi-04-1226
Fig. 7
165mm



AMi-04-1227
Fig. 5
145mm



AMi-04-1228
Fig. 5A
145mm



AMi-04-1229
Fig. 3L
210mm



AMi-04-1230
Fig. 4L
218mm




AMi-04-1231
Double Handle
No.4 + No.3
163mm



- Инструменты для амальгамы
- Instruments pour amalgame
- AmalgamisInstrumente
- Инструменти для амальгами
- Instrumentos para amalgama
- الصكوك الملعقم

AMALGAM CARRIERS



AMi-1240
Distal
Improves visibility and access to posterior regions.

AMi-1241
Double-End

AMi-1242
Mini/
Regular

AMi-1243
Large/
Regular

AMi-1244
Jumbo/
Regular

TIP SIZES FOR DOUBLE-END

AMi-1245 - 1.5mm-2.0mm AMi-1248 - 1.0mm-2mm
AMi-1246 - 2.0mm-2.8mm AMi-1249 - 2.0mm-3mm
AMi-1247 - 2.0mm-3.2mm AMi-1250 - 2.0mm-4mm

AMi1-1250 - 1.5mm-2.0mm AMi4-1250 - 1.0mm-2mm
AMi2-1250 - 2.0mm-2.8mm AMi5-1250 - 2.0mm-3mm
AMi3-1250 - 2.0mm-3.2mm AMi5-1250 - 2.0mm-4mm

AMi7-1250 - 1.5mm-2.0mm AMi10-1250 - 1.0mm-2mm
AMi8-1250 - 2.0mm-2.8mm AMi11-1250 - 2.0mm-3mm
AMi8-1250 - 2.0mm-3.2mm AMi12-1250 - 2.0mm-4mm

AMi13-1250 - 1.5mm-2.0mm AMi16-1250 - 1.0mm-2mm
AMi14-1250 - 2.0mm-2.8mm AMi17-1250 - 2.0mm-3mm
AMi15-1250 - 2.0mm-3.2mm AMi18-1250 - 2.0mm-4mm

AMi19-1250 - 1.5mm
AMi20-1250 - 2.0mm
AMi21-1250 - 2.0mm
AMi22-1250 - 3.0mm

AMi-1251
Amalgam
Well



ROOT PLANNING FILES



AMi-1252
Kramer-Nevins
Bone Chisel



AMi-1253
Kramer-Nevins
Bone Chisel



AMi-1254
Molt 18.5cm
Auto Clavable
Plastic Tips.



AMi-1255
Mead 18.5cm
Non Auto Clavable
Plastic Tips.



AMi-1256
Miller



AMi-1257
Miller

BONE FILES



AMi-1258
Miller 21



AMi-1259
Miller 1X
cross cut



AMi-1260
Miller 2X
cross cut



AMi-1261
Miller 52



AMi-1262
Wahl 1



- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Инструменти пародонтологічні
- Instrumentos para periodoncia
- الصكوك اللثة

PERIOSTEEL ELEVATORS





CHISELS



AMi-04-1232 AMi-04-1233 AMi-04-1234 AMi-04-1235



AMi-04-1236
4mm (Straight)



AMi-04-1237
4mm (Curved)



AMi-04-1238
6mm (Straight)



AMi-04-1239
6mm (Curved)

Periodontal Chisels

Bone Splitting Chisels

2.0mm

2.5mm

3.0mm

3.5mm

4.0mm

4.5mm

5.0mm



AMi-04-1240



AMi-04-1241



AMi-04-1242



AMi-04-1243



AMi-04-1244



AMi-04-1245



AMi-04-1246



AMi-04-1247
Metal tip
220 mm
Hammer With Wooden Handle



AMi-04-1248
Metal tip
220 mm



AMi-04-1249
1
For plaster



AMi-04-1250
2
For metal



AMi-04-1251
Fret saw
250 mm



AMi-04-1252
Small
Saw Blades



AMi-04-1253
Large
Saw Blades



AMi-04-1254
Ceri Saw



AMi-04-1255
Plaster Shears
200 mm



Crucible Tongs
AMi-04-1256 250 mm
AMi-04-1257 320 mm



AMi-04-1258
Ceri Saw
180 mm



AMi-04-1259
Blade
0,12x5 mm
Conf. 12 pcs



AMi-04-1260
Blade
0,25x5 mm
Conf. 12 pcs



AMi-04-1261
Blade
0,25x3 mm
Conf. 12 pcs