

HANDLES OPTIONS FOR ROOT ELEVATORS A, B, C, D, E, F, G





Dental Root Elevators

To look at one of these instruments you might mistake it for a narrow screwdriver.

That's because just like one an elevator has a handle and then a specially designed "blade" or tip portion.

How they're used.

The tip of an elevator is designed to be wedged into the ligament space between a tooth and its surrounding bone (as shown in our animation).


As the elevator is forced into and twisted around in this space, the tooth is rocked around and pressed against the walls of its socket.

This helps to both expand the shape of the socket and separate the tooth from its ligament.

As this work is continued, the tooth gradually becomes more and more mobile.

In some cases, the dentist may be able to completely remove the tooth just using their elevator.

If not, they will choose a point where they feel they have accomplished as much as they can with it and will then switch to using extraction forceps.



ROOT ELEVATOR

ASIA CARE + TECH



- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Слеватор корневий екстракційний
- Elevadores para raices
- المصاعد الجذر

Apexo Root Elevators



Apexo



AMi-150
Apexo
2.5x9mm
Straight



AMi-151
Apexo
2.7x9mm
Straight



AMi-152
Apexo
3.0x9mm
Straight



AMi-153
Apexo
2.5x9mm
Curved



AMi-154
Apexo
3.0x9mm
Curved



AMi-155
Seldin
3.0x14mm
Straight

Luxator Root Elevators



Luxator



AMi-156
Luxator
3mm
straight



AMi-157
Luxator
3mm
curved



AMi-158
Luxator
4mm
straight



AMi-159
Luxator
4mm
curved



AMi-160
Luxator
5mm
straight



AMi-161
Luxator
5mm
curved

Bein Root Elevators




Bein



AMi-162
Bein 2mm



AMi-163
Bein 3mm



AMi-164
Bein 4mm

Seldin Root Elevators



Seldin



AMi-165
Flohr 1



AMi-166
Seldin 34S



AMi-167
Seldin 34





- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber

- Elevator corneveiy ekstraktsionniy
- Elevadores para raices
- المصاعد الجذر

Apexo Root Elevators



Apexo

AMi-168
Apexo
2.5x9mm
Straight



AMi-169
Apexo
2.7x9mm
Straight



AMi-170
Apexo
3.0x9mm
Straight



AMi-171
Apexo
2.5x9mm
Curved



AMi-172
Apexo
3.0x9mm
Curved



AMi-173
Seldin
3.0x14mm
Straight



Luxator Root Elevators



Luxator

AMi-174
Luxator
3mm
straight



AMi-175
Luxator
3mm
curved



AMi-176
Luxator
4mm
straight



AMi-177
Luxator
4mm
curved



AMi-178
Luxator
5mm
straight



AMi-179
Luxator
5mm
curved




Bein Root Elevators



Bein

AMi-180
Bein 2mm



AMi-181
Bein 3mm



AMi-182
Bein 4mm



Seldin Root Elevators



Seldin

AMi-183
Flohr 1



AMi-184
Seldin 34S



AMi-185
Seldin 34




ROOT ELEVATOR

ASIA CARE + TECH



- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Єлеватор корневий екстракційний
- Elevadores para raices
- المصاعد الجذر

Apical Root Elevators




AMi-186
Fig. 1
str.

AMi-187
Fig. 4
Left

AMi-188
Fig. 5
Right

Apical Root Elevators




AMi-189
Fig. 2
str.

AMi-190
Fig. 3
Left

AMi-191
Fig. 4
Right

Coupland Root Elevators




AMi-192
Coupland No.1

AMi-193
Coupland No.2

AMi-194
Coupland No.3

Seldin Root Elevators



AMi-195
Seldin 4L
(Left)

AMi-196
Seldin 4R
(Right)

AMi-197
Seldin 1L
(Left)

AMi-198
Seldin 1R
(Right)

Seldin 13-3
Str.

AMi-199
Seldin
(Left)

AMi-200
Seldin
(Right)





- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber

- Elevator corneviu extracționiu
- Elevadores para raices
- المصاعد الجذر

Apical Root Elevators

Apical Root Elevators

AMi-201
Fig. 1
Str.

AMi-202
Fig. 4
Left

AMi-203
Fig. 5
Right

AMi-204
Fig. 2
Str.

AMi-205
Fig. 3
Left

AMi-206
Fig. 4
Right

Coupland
Root Elevators

Seldin
Root Elevators

AMi-207
Coupland No.1

AMi-208
Coupland No.2

AMi-209
Coupland No.3

AMi-210
Seldin 4L
(Left)

AMi-211
Seldin 4R
(Right)

AMi-212
Seldin 1L
(Left)

AMi-213
Seldin 1R
(Right)

Coupland

Seldin
13-3
Str.

AMi-214
Seldin
(Left)

AMi-215
Seldin
(Right)




ROOT ELEVATOR

ASIA CARE + TECH




- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Слеватор корневий екстракційний
- Elevadores para raices
- المصاعد الجذر

Apical Root Elevators



Apical Root Elevators



Apical Root Elevators



Laxutor Root Elevators





- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber


- Elevator cornevier extracționii
- Elevadores para raices
- المصاعد الجذر

Lindo Levian Root Elevators



Flohr Root Elevators

Heidbrink Root Elevators



ROOT ELEVATOR

ASIA CARE + TECH



- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Слеватор корневий екстракційний
- Elevadores para raices
- المصاعد الجذر

Friedmann Roots Elevator



AMi-251
Friedmann 4F
Str.

AMi-252
Friedmann 5F
Left

AMi-253
Friedmann 6F
Right

Heidbrink Roots Elevator



AMi-254
Heidbrink
Straight

AMi-255
Heidbrink
Left

AMi-256
Heidbrink
Right

Wisdom Roots Elevator



AMi-257
Wisdom
for 18 and 38

AMi-258
Wisdom
for 28 and 48

Bayonet Roots Elevator Set



AMi-259
Left

AMi-260
Right

AMi-261
bayonet

AMi-262
Curved

AMi-263
Straight

Luxator Roots Elevator



AMi-264
Luxator
Straight

AMi-265
Luxator
Curved





- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber

- Єлеватор кореневий екстракційний
- Elevadores para raices
- المصاعد الجذر

Roots Elevator

Luxier Roots Elevator



AMi-266 60

AMi-267 61

AMi-268 62

AMi-269 63

AMi-270 Luxier 3mm, Straight

AMi-271 Luxier 3mm, Curved

AMi-272 Luxier 5mm, Stright

AMi-273 Luxier 3mm, Angular

Luxier Roots Elevator

Claw-Lift Roots Elevator



AMi-274 Luxier Cvd.

AMi-275 Luxier Str.

AMi-276 Claw-Lift Right

AMi-277 Claw-Lift Left



ROOT ELEVATOR

ASIA CARE + TECH



- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Слеватор корневий екстракційний
- Elevadores para raices
- المصاعد الجذر

Hylin Roots Elevator

Cryer Roots Elevator



Kopp Roots Elevator



ROOT ELEVATOR

ASIA CARE + TECH




- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Слеватор корневий экстракційний
- Elevadores para raíces
- المصاعد الجذر

Medan Roots Elevator



Medan Roots Elevator



Roots Elevator





- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber

- Elevator corneviu extracționii
- Elevadores para raices
- المصاعد الجذر

Kopp Roots Elevator



AMi-295
Kopp
Right

AMi-296
Kopp
Left

Heber Roots Elevator



AMi-297
Heber
Bayonet

AMi-298
Heber
Bayonet


Bein Roots Elevator



AMi-299
Bein,
3mm

AMi-300
Bein,
4mm

Bein Roots Elevator



AMi-302
2mm
straight
serrated

AMi-303
3mm
straight

AMi-304
4mm
Curved

AMi-306
2mm
straight
serrated

AMi-307
3mm
straight

AMi-308
4mm
Curved

AMi-310
2mm
straight
serrated

AMi-311
3mm
straight

AMi-312
4mm
Curved

AMi-301
Bein
2mm, curved
serrated

AMi-305
Bein
2mm, curved
serrated

AMi-309
Bein
2mm, curved
serrated



ROOT ELEVATOR

ASIA CARE + TECH



- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Єлеватор коренвий екстракційний
- Elevadores para raíces
- المصاعد الجذر

Atraumatic Extraction Bayonet Root Elevators



AMi-27-294
3mm Luxating
Elevator Curved



AMi-28-294
3mm Luxating
Elevator Straight



AMi-29-294
5mm Luxating
Elevator Curved



AMi-30-294
5mm Luxating
Elevator Straight



AMi-31-294
2mm Grooved
Luxating Elevator



AMi-32-294
3mm Grooved
Luxating Elevator



AMi-33-294
4mm Grooved
Luxating Elevator



AMi-34-294
5mm Grooved
Luxating Elevator



ROOT ELEVATOR

ASIA CARE + TECH



- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Єлеватор коренвий екстракційний
- Elevadores para raíces
- المصاعد الجذر

Luxating Elevator With Periotome Tip



AMi-35-294
3mm Luxating Elevator / Periotome



AMi-36-294
4mm Luxating Elevator / Periotome



AMi-37-294
5mm Luxating Elevator / Periotome



AMi-38-294
Bayonet Elevator Straight

Atraumatic Extraction Bayonet Root Elevators



AMi-39-294
Bayonet Elevator Offset "In" Small



AMi-40-294
Bayonet Elevator Offset "Out" Small



AMi-41-294
Bernard Bayonet Elevator Offset Left



AMi-42-294
Bernard Bayonet Elevator Offset Right



ROOT ELEVATOR

ASIA CARE + TECH



- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Elevator corneviu extracționii
- Elevadores para raices
- المصاعد الجذر

Seldin Roots Elevator



- AMi-313** Cryer 27 (Left)
- AMi-314** Cryer 28 (Right)
- AMi-315** Seldin #1L
- AMi-316** Seldin #1R
- AMi-317** Seldin #4L
- AMi-318** Seldin #4R
- AMi-319** Seldin No.1 (Stout)
- AMi-320** Seldin No.2
- AMi-321** Kopp No.3
- AMi-322** Kopp No.4

W. Barry Roots Elevator



- AMi-323** W. Barry 1 (mm 125)
- AMi-324** W. Barry 1
- AMi-325** W. Barry 2
- AMi-326** W. Barry 3
- AMi-327** W. Barry 4

Winter Roots Elevator



- AMi-328** Winter 1L (mm 125)
- AMi-329** Winter 1L
- AMi-330** Winter 1R
- AMi-331** Winter 11L
- AMi-332** Winter 11R
- AMi-333** Winter 12L
- AMi-334** Winter 12R

Pott Roots Elevator



- AMi-335** Lecluse bayonet
- AMi-336** Pott Left
- AMi-337** Pott Right





- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber

- Єлеватор кореневий екстракційний
- Elevadores para raices
- المصاعد الجذر

Heidbrink Roots Elevator

Syndesmotomes Roots Elevator



AMi-338
Heidbrink
Straight

AMi-339
Heidbrink
Left

AMi-340
Heidbrink
Right

AMi-341
Heidbrink
Straight

AMi-342
Heidbrink
Left

AMi-343
Heidbrink
Right

AMi-344
Straight

AMi-345
Left

AMi-346
Right

Periotome Roots Elevator

Syndesmotomes



AMi-347
Periotome
2.1mm
Straight

AMi-348
Periotome
2.6mm
Straight

AMi-349
Periotome
2.1mm
Angular

AMi-350
Periotome
2mm
Curved

AMi-351
Periotome
2.0mm

AMi-352
Periotome
2.5mm

AMi-353
Periotome
3.0mm

AMi-354
Periotome
4.0mm





- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber

- Єлеватор корневий екстракційний
- Elevadores para raices
- المصاعد الجذر



FUNCTION: To loosen small root fragments from bony socket



Heidbrink Roots Elevator

Warwick-James



AMi-44-294 Heidbrink Root Apex Elevator (Left)
AMi-45-294 Heidbrink Root Apex Elevator (Right)



AMi-46-294

AMi-47-294

AMi-48-294



AMi-49-294 Fig. 2



AMi-50-294 Fig. 1



AMi-51-294 Fig. 3

AMi-43-294 Heidbrink Root Apex Elevator

ROOT ELEVATOR

ASIA CARE + TECH



- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Єлеватор кореневий екстракційний
- Elevadores para raices
- المصاعد الجذر



AMi-355
1S
(Straight 1mm)



AMi-356
2S
(Straight 2mm)



AMi-357
3S
(Straight 3mm)



AMi-358
5S
(Straight 5mm)



AMi-359
3C
(Curved 3mm)



AMi-360
5C
(Curved 5mm)



AMi-361
3CA
(Curved Angle 3mm)



AMi-362
3IC
(Inverted Curved 3mm)



AMi-363
DE 3
(Dual Edge 3mm / 1.5mm)



AMi-364
DE 5
(Dual Edge 5mm / 3mm)





- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber

- Єлеватор корневий екстракційний
- Elevadores para raices
- المصاعد الجذر



AMi-52-294
(Straight 1.5mm)



AMi-53-294
(Curved 1.5mm)



AMi-54-294
(Straight 2mm)



AMi-55-294
(Curved 2mm)



AMi-56-294
(Straight 3mm)



AMi-57-294
(Curved 3mm)



AMi-58-294
(Straight 4mm)



AMi-59-294
(Curved 4mm)



AMi-60-294
(Straight 5mm)



AMi-61-294
(Curved 5mm)

AMi-62-294
PERIOTOME LUXATING ELEVATORS SET OF 04 PCS - GOLD COATED TIPS



AMi-63-294
Periotome Gold Tip - L3S



AMi-64-294
Periotome Gold Tip - De3



AMi-65-294
Periotome Gold Tip - L2S




AMi-66-294
Periotome Gold Tip - L3CA




- Элеватор корневой экстракционный
 - **Elévateurs à racines**
 - Wurzelheber
- PDL ELEVATORS**

- Єлеватор кореневий екстракційний
- Elevadores para raices
- المصاعد الجذر



AMi-365
PDL Elevators
Set of 7 Pcs



AMi-366
PDL Elevators
Set of 7 Pcs

ROOT ELEVATOR

ASIA CARE + TECH



- Элеватор корневой экстракционный
- Elévateurs à racines
- Wurzelheber
- Єлеватор корневий екстракційний
- Elevadores para raices
- المصاعد الجذر

PDL ELEVATORS

Delicate Arrow tip

Luxivators straight
atraumatic luxation of the tooth

Luxivators curved
atraumatic luxation of the tooth



4,0mm straight
4,0mm curved

2,5mm straight
3,5mm straight
4,5mm straight

2,5mm curved
3,5mm curved
4,5mm curved

AMi-367
PDL Elevators
Set of 8 Pcs

